
# Adult Cardiac Arrest Circular Algorithm Assessment and Treatment

# Circular Algorithm





### **CPR**

#### **CPR Quality:**

- Push hard (at least 2 inches [5 cm]) and fast (100-120/min) and allow complete chest recoil.
- Minimize interruptions in compressions.
- Avoid excessive ventilation.
- Change compressor every 2 minutes, or sooner if fatigued.
- If no advanced airway, 30:2 compression-ventilation ratio, or 1 breath every 6 seconds.
- Quantitative waveform capnography
- If PETCO2 is low or decreasing, reassess CPR quality.

#### **Shock Energy for Defibrillation:**

- Biphasic: Manufacturer recommendation (eg, initial dose of 120-200 J); if unknown, use maximum available. Second and subsequent doses should be equivalent, and higher doses may be considered.
- · Monophasic: 360 J

#### **Drug Therapy:**

- Epinephrine IV/IO dose: 1 mg every 3-5 minutes
- Amiodarone IV/IO dose: First dose: 300 mg bolus. Second dose: 150 mg.
- Lidocaine IV/IO dose: First dose: 1-1.5 mg/kg. Second dose: 0.5-0.75 mg/kg.

# **Advanced Airway:**

- Endotracheal intubation or supraglottic advanced airway
- Waveform capnography or capnometry to confirm and monitor ET tube placement
- Once advanced airway in place, give 1 breath every 6 seconds (10 breaths/min) with continuous chest compressions

#### **Return of Spontaneous Circulation (ROSC):**

- Pulse and blood pressure
- Abrupt sustained increase in PETCO2 (typically ≥40 mm Hg)
- Spontaneous arterial pressure waves with intra-arterial monitoring

## **Reversible Causes:**

- Hypovolemia
- Hypoxia
- Hydrogen ion (acidosis)
- Hypo-/hyperkalemia
- Hypothermia
- Tension pneumothorax
- Tamponade, cardiac
- Toxins
- Thrombosis, pulmonary
- Thrombosis, coronary